

## Introducing the First and Only Oxygen and Air Flowmeters with Antimicrobial Product Protection

The antimicrobial-treated rotary knob on these flowmeters will continue to inhibit the growth of odor causing bacteria, mold, mildew, and fungi on the surface of the knobs, and maintain antimicrobial efficacy for the lifespan of the product.

- Built-in product protection prevents growth of microbes on the surface of the product
- Broad Regulatory Approvals including EPA, FDA, NSF, Health Canada, European BDP, Swiss & EFSA
- Gentec flowmeters have a 0-15 LPM flow range, with a minimum of 65 LPM flush flow rate
- Antimicrobial based on non-volatile, naturally-occurring silver
- Antimicrobial long-lasting product protection, lasts the lifetime of the knob
- Environmentally-friendly technology (green product)
- Cradle-to-cradle silver certified
- Available in either nickel plated solid brass or aluminum
- Durable polycarbonate hood and metering tube
- Stainless steel ball float
- Flowmeters are all calibrated at 50 psig and 70° F
- Flowmeters are all back-pressured compensated
- Color-coded for air and oxygen
- 5-year warranty (see warranty statement for details)



Nature's antimicrobial

Agion® is a registered trademark of Agion® Technologies.



FM197B-15L-CH-P  
Oxygen Flowmeter w/ Power Take-off & Chemetron Adapter

FM187B-15L-00

## FLOWMETER ORDERING INFORMATION

Flowmeter Type  Adapter	Air Flowmeter 0-15 LPM	Air Flowmeter w/ Power Take-off	Air Dual Flowmeter 0-15 LPM	Oxygen Flowmeter 0-15 LPM	Oxygen w/ Power Take-off	Oxygen Dual Flowmeter 0-15 LPM
Standard Inlet 1/8" NPT Female	FM187B-15L-00-E	FM187B-15L-00-EP	FMD187B-15L-00-E	FM197B-15L-00-E	FM197B-15L-00-EP	FMD197B-15L-00-E
OHMEDA Adapter	FM187B-15L-OH-E	FM187B-15L-OH-EP	FMD187B-15L-OH-E	FM197B-15L-OH-E	FM197B-15L-OH-EP	FMD197B-15L-OH-E
CHEMETRON Adapter	FM187B-15L-CH-E	FM187B-15L-CH-EP	FMD187B-15L-CH-E	FM197B-15L-CH-E	FM197B-15L-CH-EP	FMD197B-15L-CH-E
DISS Hand-Tight Nut & Nipple Adapter	FM187B-15L-DH-E	FM187B-15L-DH-EP	FMD187B-15L-DH-E	FM197B-15L-DH-E	FM197B-15L-DH-EP	FMD197B-15L-DH-E
Diss Nut & Nipple Adapter	FM187B-15L-DS-E	FM187B-15L-DS-EP	FMD187B-15L-DS-E	FM197B-15L-DS-E	FM197B-15L-DS-EP	FMD197B-15L-DS-E
Puritan-Bennett Adapter	FM187B-15L-PB-E	FM187B-15L-PB-EP	FMD187B-15L-PB-E	FM197B-15L-PB-E	FM197B-15L-PB-EP	FMD197B-15L-PB-E
Oxequip Medstar Adapter	FM187B-15L-OE-E	FM187B-15L-OE-EP	FMD187B-15L-OE-E	FM197B-15L-OE-E	FM197B-15L-OE-EP	FMD197B-15L-OE-E
Oxequip Twist Adapter	FM187B-15L-OQ-E	FM187B-15L-OQ-EP	FMD187B-15L-OQ-E	FM197B-15L-OQ-E	FM197B-15L-OQ-EP	FMD197B-15L-OQ-E

**OPTIONS:**

- For aluminum flowmeter, change the B to an A (example: FM197B-15L-OH to FM197A-15L-OH)
- For dual flowmeter with 2 power take-offs, add -2P (example: FMD197B-15L-OH-2P)
- For ISO flowmeters add an E (example: FM197A-15L-OHE)

## SPECIFICATIONS

<b>Weight (no adapter)</b>	Aluminum Body: 4.24 oz/122 gm Brass Body: 6.64 oz/188 gm
<b>Increment:</b>	0-6 LPM in ½ LPM increment, 6-15 LPM in 1 LPM Increment
<b>Flush Flow Rate:</b>	65 LPM (depending on supply flow)
<b>Max Inlet Pressure:</b>	100 psig
<b>PTO Flow Rate:</b>	Over 100 LPM (depending on supply flow)