



## **DIGITAL SUCTION REGULATORS**

Gentec is proud to announce the release of the all-new Digital Suction Regulator. While preserving the look-and-feel and usability of our original 2 ½" gauge analog version, we provide the added convenience of an easy-to-read digital display (including a digital sweeping needle.) Truly the best of both worlds!.

## Agion®\* Antimicrobial Product Protection

The built-in long-lasting product protection prevents growth of microbes on the surface of the knobs. The antimicrobial product protection is based on non-volatile, naturally occurring silver.





- Lightweight with durable ABS plastic casing
- Impact-resistant internal components
- Corrosion-free components and can be thoroughly flushed
- On/off cycle settings can be adjusted without opening the case
- Accuracy: +/-1% for digital; +/- 2% for analog
- Intermittent timing cycles: 16 seconds on, 8 seconds off. (timing is adjustable)
- Suction traps and a variety of fittings & adapters are available
- Color-coding available
- 10 year warranty




Agion's silver antimicrobial has been certified for its material content, recyclability, and manufacturing characteristics. Cradle To Cradle Certified<sup>CM</sup> products meet established standards for human health, environmental health and recyclability. Cradle To Cradle Certification mark of MBDC. **Digital Suction Regulators** 



## **ORDERING MATRIX**

Instructions: Select the suction regulator, adapter, fitting, and color (if needed), then write the model number below.



Medical Gas Resources 5592 Buckingham Dr., Huntington Beach, CA 92649 T. 888.743.8298 F. 714.897.7426 www.medicalgasresources.com